

Classification of gastrointestinal lesions for resection

Project group 10

Athena Taymourtash


Bodnár Dániel

Alice Santilli

Objective – problem specification


❖ Colorectal cancer and the importance of the accuracy and the promptness of the diagnosis

- Cause: Old age, lifestyle factors, genetic disorders
- Diagnostic method: Tissue biopsy and colonoscopy



❖ Virtual Colonoscopy

- tell surgeons in real time whether or not the lesion needs to be resected
- discriminate the severity of individual lesions in patients with many polyps in real time



Examples of the available images

- videos, 25 frames/second, 768x576 pixels per frame
- each video length is 5 - 30 seconds
- 76 colonoscopy videos, both in NBI (narrow band imaging) and WL (white light)
- 15 serrated adenomas, 21 hyperplastic lesions and 40 adenomas
- Labels :

Adenoma + serrated = resection

hyperplastic = no resection

	Human Expert	Human Beginner
Acc.	76.54%	72.22%
Sen.	64.98%	54.50%
Spec.	81.35%	77.30%



















































Model Architecture

- By considering a video as a series of frames:
 1. Loop over all extracted frames in the video file
 2. For each frame, pass the frame through the CNN
 3. Classify each frame *individually* and *independently* of each other
 4. Choose the label with the largest corresponding probability
- For test data:
 - Computing the average of the predictions for all frames of a video and choosing the label with the largest corresponding probability

Data separation:
















Training

(D:) > Gastro-Intestinal-Lesions > data > train

 adenoma_01_WL	 hyperplasic_08_WL	 adenoma_28_WL
 adenoma_02_WL	 hyperplasic_09_WL	 adenoma_29_WL
 adenoma_03_WL	 hyperplasic_13_WL	 adenoma_34_WL
 adenoma_04_WL	 hyperplasic_14_WL	 adenoma_37_WL
 adenoma_05_WL	 hyperplasic_15_WL	 hyperplasic_01_WL
 adenoma_06_WL	 hyperplasic_16_WL	 hyperplasic_02_WL
 adenoma_07_WL	 hyperplasic_18_WL	 hyperplasic_03_WL
 adenoma_08_WL	 hyperplasic_19_WL	 hyperplasic_04_WL
 adenoma_09_WL	 serrated_01_WL	 hyperplasic_06_WL
 adenoma_13_WL	 serrated_02_WL	 hyperplasic_07_WL
 adenoma_14_WL	 serrated_03_WL	
 adenoma_15_WL	 serrated_04_WL	
 adenoma_16_WL	 serrated_05_WL	
 adenoma_17_WL	 serrated_06_WL	
 adenoma_18_WL	 serrated_07_WL	
 adenoma_19_WL	 serrated_08_WL	
 adenoma_20_WL	 serrated_13_WL	
 adenoma_22_WL	 serrated_14_WL	
 adenoma_24_WL	 serrated_15_WL	


Validation

(D:) > Gastro-Intestinal-Lesions > data > val


-  adenoma_10_WL
-  adenoma_11_WL
-  adenoma_12_WL
-  adenoma_25_WL
-  adenoma_30_WL
-  adenoma_31_WL
-  adenoma_39_WL
-  adenoma_40_WL
-  hyperplasic_10_WL
-  hyperplasic_11_WL
-  hyperplasic_12_WL
-  hyperplasic_20_WL
-  hyperplasic_21_WL
-  serrated_10_WL
-  serrated_12_WL

Test

(D:) > Gastro-Intestinal-Lesions > data > test



adenoma_21_WL



adenoma_23_WL




adenoma_26_WL




hyperplasic_05_WL



hyperplasic_17_WL



serrated_09_WL



serrated_11_WL

Preprocessing

Colonoscopy_Train.py

```
train_data, train_labels = DataLoader("train", datapath=args.dataset, step=10)
val_data, val_labels = DataLoader("val", datapath=args.dataset, step=10)
```

```
def DataLoader(what, datapath, step):
```

```
    STEP = step
    frames = []
    labels = []
    f_counter = 0
```

```
    if what == "train":
        print("loading training data.....")
        directory = datapath + 'train/'
```

```
    if what == "val":
        print("loading validation data.....")
        directory = datapath + 'val/'
```

```
    for filename in os.listdir(directory):
```

```
        video = cv2.VideoCapture(directory + filename)
        name = filename.split("_")[0] #get the class of the video
        if name=="adenoma" or name=="serrated":
            label = [0,1]
        elif name=='hyperplasic':
            label = [1,0] #one hot encoded labels
        else :
            print("issue with filename " + filename)
```

```
    while True:
        has_frame, frame = video.read()
```

```
        if not has_frame:
            break
```

```
        if f_counter % STEP == 0:
            frames.append(frame)
            labels.append(label)
```


```
        f_counter += 1
```

Training Model

- Transfer learning of ResNet50
- Froze all the base layers
- Fine tuned training to our dataset

Retrain ResNet50

ResNet50 Diagram



Training Performance

- Short training time due to transfer learning
- 80% validation accuracy
- 97% training accuracy



Results

Predictions were made for the entire video

Test on 7 full videos :

5 needing resection

(3 adenoma, 2 serrated)

2 no resection

(hyperplastic polyp)

Type	True Class	Predicted Class	Confidence
Adenoma	1	1	100%
Adenoma	1	1	98%
Adenoma	1	1	99%
Serrated	1	1	97%
Serrated	1	1	51%
Hyperplastic	0	0	63%
Hyperplastic	0	0	56%

Future work

- Give surgeons information about which type of polyp they are viewing
 - This may help them make better clinical decisions
- Try different architectures or help extract more features from the images to aid the classification



Hyperplastic



Adenoma



Serrated

Conclusion

Colon Polyps and lesions can be successfully identified from colonoscopy videos.

This will give surgeons real time pathologic information about the polyps and make sure all important lesions are being resected.

Thank you